

LOX FILL HARNESS KIT

P/N 50C-0097-10
P/N 50C-0097-15



LOX Fill Harness Kit connects a Liquid Oxygen Cylinder (VGL) with a LOX converter or system for filling. The kit includes a relief valve, manual vent valve, CRU-59 LOX filler valve and hose assembly in either a 10' or 15' length.

DIMENSIONS (W" x L" x H")	23 x 18.5 x 7
WEIGHT	15 lbs
RELIEF VALVE SETTING	400 PSIG
SUPPLY CONNECTION	CRU-59/E Fill Valve

BACKPACK MEDICAL OXYGEN SYSTEM FILLING STATION HOSE KIT

P/N 50C-0116-1

BMOS-FS Hose Kit connects the MMOS to the BMOS-FS for refilling. The kit includes a relief valve, burst disc, CRU-50 LOX filler valve and CRU-59 LOX filler valve attached to a 4' cryogenic transfer hose. Items are packaged in a protective storage case.

DIMENSIONS (W" x L" x H")	18.5 x 14 x 7
WEIGHT	13 lbs
RELIEF VALVE SETTING	75 PSIG
SAFETY BURST DISC PRESSURE	150 ± 10 PSIG
OUTLET CONNECTION	CRU-59/E Fill Valve on a 4' Hose
SUPPLY CONNECTION	CRU-50/A Fill Fitting

